

For Research Use Only

BTBD6 Polyclonal antibody

Catalog Number: 13776-1-AP



Basic Information

Catalog Number:

13776-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4766

GenBank Accession Number:

BC042525

GeneID (NCBI):

90135

UNIPROT ID:

Q96KE9

Full Name:

BTB (POZ) domain containing 6

Calculated MW:

410 aa, 46 kDa

Observed MW:

53-58 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HepG2 cells,

IP : Y79 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 13776-1-AP (BTBD6 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.

IP result of anti-BTBD6 (IP:13776-1-AP, 4ug; Detection:13776-1-AP 1:300) with Y79 cells lysate 800ug.