

For Research Use Only

WBSCR16 Polyclonal antibody

Catalog Number: 13796-1-AP

Featured Product

4 Publications



Basic Information

Catalog Number:

13796-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4619

GenBank Accession Number:

BC032712

GeneID (NCBI):

81554

UNIPROT ID:

Q96I51

Full Name:

Williams-Beuren syndrome
chromosome region 16

Calculated MW:

464 aa, 50 kDa

Observed MW:

50 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : HeLa cells, Jurkat cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Aurelio Reyes	32735630	PLoS Genet	WB
Hana Antonicka	27974379	EMBO Rep	WB
Nicholas J Kramer	37770567	Nat Cell Biol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 13796-1-AP (WBSCR16 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.