

For Research Use Only

SLITRK4 Polyclonal antibody

Catalog Number: 13835-1-AP **1 Publications**



Basic Information

Catalog Number: 13835-1-AP	GenBank Accession Number: BC040986	Purification Method: Antigen affinity purification
Size: 193 µg/ml	GeneID (NCBI): 139065	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	UNIPROT ID: Q8IW52	
Isotype: IgG	Full Name: SLIT and NTRK-like family, member 4	
Immunogen Catalog Number: AG4781	Calculated MW: 837 aa, 94 kDa	
	Observed MW: 96 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse thymus tissue,
Cited Applications: WB	
Species Specificity: human, mouse, rat	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yao-Qi Zhou	37012514	Cell Oncol (Dordr)	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 13835-1-AP (SLITRK4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.