

For Research Use Only

BRUNOL6 Polyclonal antibody

Catalog Number: 13866-1-AP



Basic Information

Catalog Number: 13866-1-AP	GenBank Accession Number: BC030835	Purification Method: Antigen affinity purification
Size: 310 µg/ml	GeneID (NCBI): 60677	
Source: Rabbit	UNIPROT ID: Q96J87	
Isotype: IgG	Full Name: bruno-like 6, RNA binding protein (Drosophila)	
Immunogen Catalog Number: AG4848	Calculated MW: 481 aa, 50 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data