

For Research Use Only

PCDHB12 Polyclonal antibody

Catalog Number: 14020-1-AP



Basic Information

Catalog Number:

14020-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5135

GenBank Accession Number:

BC045637

GeneID (NCBI):

56124

UNIPROT ID:

Q9Y5F1

Full Name:

protocadherin beta 12

Calculated MW:

87 kDa

Observed MW:

87 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse brain tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 14020-1-AP (PCDHB12 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.