

For Research Use Only

MORC1 Polyclonal antibody

Catalog Number: 14080-1-AP

Featured Product

1 Publications



Basic Information

Catalog Number:

14080-1-AP

Size:

287 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5205

GenBank Accession Number:

BC050307

GeneID (NCBI):

27136

UNIPROT ID:

Q86VD1

Full Name:

MORC family CW-type zinc finger 1

Calculated MW:

112 aa, 8 kDa

Observed MW:

108 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2400

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

rat, mouse

Positive Controls:

WB : mouse testis tissue, human testis tissue

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|--------------------|-----------|---------------|-------------|
| Annakarina Mundorf | 34331083 | Exp Brain Res | WB |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 14080-1-AP (MORC1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.