

For Research Use Only

HTRA4 Polyclonal antibody

Catalog Number: 14095-1-AP



Basic Information

Catalog Number:

14095-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5227

GenBank Accession Number:

BC057765

GeneID (NCBI):

203100

UNIPROT ID:

P83105

Full Name:

HtrA serine peptidase 4

Calculated MW:

476 aa, 51 kDa

Observed MW:

51 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : human placenta tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 14095-1-AP (HTRA4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.