

For Research Use Only

# RAP1GAP Polyclonal antibody

Catalog Number: 14229-1-AP



## Basic Information

Catalog Number: 14229-1-AP	GenBank Accession Number: BC054490	Purification Method: Antigen affinity purification
Size: 273 µg/ml	GeneID (NCBI): 5909	Recommended Dilutions: IHC 1:20-1:200
Source: Rabbit	UNIPROT ID: P47736	
Isotype: IgG	Full Name: RAP1 GTPase activating protein	
Immunogen Catalog Number: AG5469	Calculated MW: 73 kDa	

## Applications

Tested Applications:  
IHC, ELISA

Species Specificity:  
human, mouse, rat

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

Positive Controls:

IHC : human pancreas tissue,

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human pancreas tissue slide using 14229-1-AP (RAP1GAP Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human pancreas tissue slide using 14229-1-AP (RAP1GAP Antibody) at dilution of 1:50 (under 40x lens).