

For Research Use Only

NFAM1 Polyclonal antibody

Catalog Number: 14310-1-AP



Basic Information

Catalog Number: 14310-1-AP	GenBank Accession Number: BC038241	Purification Method: Antigen affinity purification
Size: 175 µg/ml	GeneID (NCBI): 150372	Recommended Dilutions: IHC 1:20-1:200
Source: Rabbit	UNIPROT ID: Q8NET5	
Isotype: IgG	Full Name: NFAT activating protein with ITAM motif 1	
Immunogen Catalog Number: AG5505	Calculated MW: 30 kDa	

Applications

Tested Applications:
IHC, ELISA

Species Specificity:
human

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IHC: human tonsillitis tissue, human lung cancer tissue, human liver cancer tissue

Background Information

NFAM1, also named as CNAIP, is an immunoreceptor tyrosine-based activation motif-bearing molecule that regulates B cell development and signaling.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 14310-1-AP (NFAM1 Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 14310-1-AP (NFAM1 Antibody) at dilution of 1:50 (under 40x lens).