

For Research Use Only

KAT2A/GCN5 Polyclonal antibody

Catalog Number: 14983-1-AP

Featured Product

5 Publications



Basic Information

Catalog Number:

14983-1-AP

Size:

210 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6940

GenBank Accession Number:

BC032743

GeneID (NCBI):

2648

UNIPROT ID:

Q92830

Full Name:

K(lysine) acetyltransferase 2A

Calculated MW:

94 kDa

Observed MW:

94 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human, mouse

Positive Controls:

WB: MCF7 cells, mouse ovary tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Bo Wang	31115516	Mol Med Rep	WB
Bowen Liu	29471853	J Hematol Oncol	WB
Leilei Li	25686500	Biochem Biophys Res Commun	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF7 cells were subjected to SDS PAGE followed by western blot with 14983-1-AP (GCN5 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.