

For Research Use Only

UNKL Polyclonal antibody

Catalog Number: 15000-1-AP **2 Publications**



Basic Information

Catalog Number:

15000-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6990

GenBank Accession Number:

BC000150

GeneID (NCBI):

64718

UNIPROT ID:

Q9H9P5

Full Name:

unkempt homolog (Drosophila)-like

Calculated MW:

25 kDa

Observed MW:

32 kDa, 45 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : HeLa cells,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Shiyong Zhao	26844847	Cancer Biother Radiopharm	WB
Wei Peng	37397926	Toxicol Res (Camb)	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 15000-1-AP (UNKL antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.