

For Research Use Only

CLCF1 Polyclonal antibody

Catalog Number: 15172-1-AP



Basic Information

| | | |
|--|---|--|
| Catalog Number: 15172-1-AP | GenBank Accession Number: BC012939 | Purification Method: Antigen affinity purification |
| Size: 110 µg/ml | GeneID (NCBI): 23529 | |
| Source: Rabbit | UNIPROT ID: Q9UBD9 | |
| Isotype: IgG | Full Name: cardiotrophin-like cytokine factor 1 | |
| Immunogen Catalog Number: AG3204 | Calculated MW: 25 kDa | |

Applications

Tested Applications:
ELISA

Species Specificity:
human, mouse, rat

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data