

For Research Use Only

TRADD Polyclonal antibody, PBS Only

Catalog Number: 15468-1-PBS

Featured Product



Basic Information

Catalog Number:

15468-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7768

GenBank Accession Number:

BC004491

GeneID (NCBI):

8717

UNIPROT ID:

Q15628

Full Name:

TNFRSF1A-associated via death domain

Calculated MW:

34 kDa

Observed MW:

34 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IP, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 15468-1-AP (TRADD antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 15468-1-PBS in a different storage buffer formulation.



IP result of anti-TRADD (IP:15468-1-AP, 4ug; Detection:15468-1-AP 1:500) with Jurkat cells lysate 4000 ug. This data was developed using the same antibody clone with 15468-1-PBS in a different storage buffer formulation.



WB result of TRADD antibody (15468-1-AP; 1:2000; incubated at room temperature for 1.5 hours) with sh-Control and sh-TRADD transfected A549 cells. This data was developed using the same antibody clone with 15468-1-PBS in a different storage buffer formulation.