

For Research Use Only

# PITPNC1 Polyclonal antibody

Catalog Number: 15715-1-AP



## Basic Information

Catalog Number:

15715-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8347

GenBank Accession Number:

BC007905

GeneID (NCBI):

26207

UNIPROT ID:

Q9UKF7

Full Name:

phosphatidylinositol transfer protein,  
cytoplasmic 1

Calculated MW:

332aa,38 kDa; 332aa,38 kDa

Observed MW:

38 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP 0.5-4.0 µg for 1.0-3.0 mg of total  
protein lysate

## Applications

Tested Applications:

IP, ELISA

Species Specificity:

human, mouse

Positive Controls:

IP: mouse brain tissue,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



IP result of anti-PITPNC1 (IP:15715-1-AP, 4ug; Detection:15715-1-AP 1:500) with mouse brain tissue lysate 5200ug.