

For Research Use Only

HLA-DRB4 Polyclonal antibody

Catalog Number: 15733-1-AP



Basic Information

Catalog Number:

15733-1-AP

Size:

260 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8492

GenBank Accession Number:

BC005312

GeneID (NCBI):

3126

UNIPROT ID:

P13762

Full Name:

major histocompatibility complex,
class II, DR beta 4

Calculated MW:

266 aa, 30 kDa

Observed MW:

35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB: Jurkat cells, human plasma tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 15733-1-AP (HLA-DRB4 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.