

For Research Use Only

UBE1L Polyclonal antibody

Catalog Number: 15818-1-AP **2 Publications**



Basic Information

Catalog Number: 15818-1-AP	GenBank Accession Number: BC006378	Purification Method: Antigen affinity purification
Size: 400 µg/ml	GeneID (NCBI): 7318	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	UNIPROT ID: P41226	
Isotype: IgG	Full Name: ubiquitin-like modifier activating enzyme 7	
Immunogen Catalog Number: AG8692	Calculated MW: 1012 aa, 112 kDa	
	Observed MW: 112 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse liver tissue,
Cited Applications: WB	
Species Specificity: human, mouse	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jie Jin	34257278	Cell Death Dis	WB
Mei-Juan Zhou	28035359	Int J Mol Med	WB


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 15818-1-AP (UBE1L antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.