

For Research Use Only

UBE1L Polyclonal antibody

Catalog Number: 15818-1-AP

2 Publications



Basic Information

Catalog Number:

15818-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8692

GenBank Accession Number:

BC006378

GeneID (NCBI):

7318

UNIPROT ID:

P41226

Full Name:

ubiquitin-like modifier activating enzyme 7

Calculated MW:

1012 aa, 112 kDa

Observed MW:

112 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : mouse liver tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Jie Jin	34257278	Cell Death Dis	WB
Mei-Juan Zhou	28035359	Int J Mol Med	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse liver tissue were subjected to SDS PAGE followed by western blot with 15818-1-AP (UBE1L antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.