

For Research Use Only

MED4 Polyclonal antibody

Catalog Number: 15951-1-AP



Basic Information

Catalog Number:

15951-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8629

GenBank Accession Number:

BC005189

GeneID (NCBI):

29079

UNIPROT ID:

Q9NPJ6

Full Name:

mediator complex subunit 4

Calculated MW:

270 aa, 30 kDa

Observed MW:

30-35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : Jurkat cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 15951-1-AP (MED4 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.