

For Research Use Only

RDH13 Polyclonal antibody

Catalog Number: 16067-1-AP

1 Publications



Basic Information

Catalog Number:

16067-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8944

GenBank Accession Number:

BC009881

GeneID (NCBI):

112724

UNIPROT ID:

Q8NBN7

Full Name:

retinol dehydrogenase 13 (all-trans/9-cis)

Calculated MW:

331 aa, 36 kDa

Observed MW:

36 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : Y79 cells, human heart tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Ying Chen	39075504	J Exp Clin Cancer Res	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Y79 cells were subjected to SDS PAGE followed by western blot with 16067-1-AP (RDH13 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.