

For Research Use Only

# UBE2Q1 Polyclonal antibody

Catalog Number: 16072-1-AP



## Basic Information

Catalog Number:

16072-1-AP

Concentration:

650 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9028

GenBank Accession Number:

BC009286

GeneID (NCBI):

55585

UNIPROT ID:

Q7Z7E8

Full Name:

ubiquitin-conjugating enzyme E2Q  
family member 1

Calculated MW:

422 aa, 46 kDa

Observed MW:

46 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP 0.5-4.0 ug for 1.0-3.0 mg of total  
protein lysate

## Applications

Tested Applications:

IP, ELISA

Species Specificity:

human

Positive Controls:

IP : HEK-293 cells,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



IP result of anti-UBE2Q1 (IP:16072-1-AP, 4ug;  
Detection:16072-1-AP 1:600) with HEK-293 cells  
lysate 1440 ug.