

For Research Use Only

LPCAT1 Polyclonal antibody

Catalog Number: 16112-1-AP

Featured Product

35 Publications



Basic Information

Catalog Number:

16112-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9035

GenBank Accession Number:

BC020166

GeneID (NCBI):

79888

UNIPROT ID:

Q8NF37

Full Name:

lysophosphatidylcholine acyltransferase 1

Calculated MW:

534 aa, 59 kDa

Observed MW:

59 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

IHC 1:50-1:500

IF/ICC 1:50-1:500

Applications

Tested Applications:

IF, IHC, IP, WB, ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, rat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB: A431 cells, human lung tissue, mouse spleen tissue, mouse brain tissue, rat brain tissue, A549 cells, mouse lung tissue, rat lung tissue

IP: mouse brain tissue,

IHC: human breast cancer tissue, mouse spleen tissue, mouse lung tissue, human colon cancer tissue

IF/ICC: MCF-7 cells,

Background Information

LPCAT1, also named as AYL2, PFAAP3 and LysoPAFAT, belongs to the 1-acyl-sn-glycerol-3-phosphate acyltransferase family. It is a key enzyme for remodeling phospholipids, including phosphatidylcholine. The expression level of LPCAT1 is able to differentiate prostate cancer from noncancerous prostatic changes, and correlates to the tumor grade of prostate cancer. LPCAT1 possesses both acyltransferase and acetyltransferase activities. It mediates the conversion of 1-acyl-sn-glycero-3-phosphocholine (LPC) into phosphatidylcholine (PC).

Notable Publications

Author	Pubmed ID	Journal	Application
Patrick Lebok	31533087	Aging (Albany NY)	IHC
Xuedan Deng	36099794	Biomed Pharmacother	IHC, WB
Martin P Helley	26434622	Neuroscience	WB, IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 16112-1-AP (LPCAT1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.




Various lysates were subjected to SDS PAGE followed by western blot with 16112-1-AP (LPCAT1 antibody) at dilution of 1:15000 incubated at room temperature for 1.5 hours.




Immunohistochemical analysis of paraffin-embedded human breast cancer tissue slide using 16112-1-AP (LPCAT1 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).




Immunohistochemical analysis of paraffin-embedded mouse spleen tissue slide using 16112-1-AP (LPCAT1 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).




Immunohistochemical analysis of paraffin-embedded mouse spleen tissue slide using 16112-1-AP (LPCAT1 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).




Immunohistochemical analysis of paraffin-embedded mouse lung tissue slide using 16112-1-AP (LPCAT1 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).





Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 16112-1-AP (LPCAT1 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 16112-1-AP (LPCAT1 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).




Immunohistochemical analysis of paraffin-embedded mouse lung tissue slide using 16112-1-AP (LPCAT1 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



IP result of anti-LPCAT1 (IP:16112-1-AP, 4ug; Detection:16112-1-AP 1:4000) with mouse brain tissue lysate 1280 ug.

Immunohistochemical analysis of paraffin-embedded human breast cancer tissue slide using 16112-1-AP (LPCAT1 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).

IP result of anti-LPCAT1 (IP:16112-1-AP, 3ug; Detection:16112-1-AP 1:4000) with mouse brain tissue lysate 3600ug.



Immunofluorescent analysis of (-20°C Ethanol) fixed MCF-7 cells using 16112-1-AP (LPCAT1 antibody) at dilution of 1:50 and Coralite488-Conjugated Goat Anti-Rabbit IgG(H+L).