

For Research Use Only

GAL3ST4 Polyclonal antibody

Catalog Number: 16240-1-AP



Basic Information

Catalog Number:

16240-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9063

GenBank Accession Number:

BC012976

GeneID (NCBI):

79690

UNIPROT ID:

Q96RP7

Full Name:

galactose-3-O-sulfotransferase 4

Calculated MW:

486 aa, 54 kDa

Observed MW:

55-60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : human placenta tissue, A2780 cells, mouse thymus tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human placenta tissue were subjected to SDS PAGE followed by western blot with 16240-1-AP (GAL3ST4 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.