

For Research Use Only

RPP21 Polyclonal antibody

Catalog Number: 16377-1-AP



Basic Information

Catalog Number: 16377-1-AP	GenBank Accession Number: BC011730	Purification Method: Antigen affinity purification
Size: 350 µg/ml	GeneID (NCBI): 79897	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q9H633	
Isotype: IgG	Full Name: ribonuclease P/MRP 21kDa subunit	
Immunogen Catalog Number: AG9353	Calculated MW: 154 aa, 18 kDa	
	Observed MW: 18 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : COLO 320 cells,
Species Specificity: human, mouse, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



COLO 320 cells were subjected to SDS PAGE followed by western blot with 16377-1-AP (RPP21 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.