

For Research Use Only

RBM38 Polyclonal antibody

Catalog Number: 16435-1-AP



Basic Information

Catalog Number:

16435-1-AP

Size:

350 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9748

GenBank Accession Number:

BC018711

GeneID (NCBI):

55544

UNIPROT ID:

Q9HOZ9

Full Name:

RNA binding motif protein 38

Calculated MW:

216 aa, 23 kDa

Observed MW:

25-31 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, rat

Positive Controls:

WB : HL-60 cells, MCF-7 cells, SMMC-7721 cells, rat heart tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 16435-1-AP (RBM38 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.