

For Research Use Only

PNMA2 Polyclonal antibody

Catalog Number: 16445-1-AP **2 Publications**



Basic Information

Catalog Number: 16445-1-AP	GenBank Accession Number: BC062301	Purification Method: Antigen affinity purification
Size: 550 µg/ml	GeneID (NCBI): 10687	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	UNIPROT ID: Q9UL42	
Isotype: IgG	Full Name: paraneoplastic antigen MA2	
Immunogen Catalog Number: AG10289	Calculated MW: 364 aa, 42 kDa	
	Observed MW: 40 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Jurkat cells, mouse brain tissue
Cited Applications: WB	
Species Specificity: human, mouse, rat	
Cited Species: mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Junjie Xu	38301645	Cell	WB
Junjie Xu	36798413	bioRxiv	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 16445-1-AP (PNMA2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.