

For Research Use Only

BEX1 Polyclonal antibody

Catalog Number: 16491-1-AP



Basic Information

Catalog Number:

16491-1-AP

Size:

110 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9616

GenBank Accession Number:

BC018615

GeneID (NCBI):

55859

UNIPROT ID:

Q9HBH7

Full Name:

brain expressed, X-linked 1

Calculated MW:

125 aa, 15 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse

Background Information

BEX1 is a signaling adapter molecule involved in p75NTR/NGFR signaling. It plays a role in cell cycle progression and neuronal differentiation. BEX1 inhibits neuronal differentiation in response to nerve growth factor (NGF).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

