

For Research Use Only

# IL-4 Polyclonal antibody

Catalog Number: 16545-1-AP **4 Publications**



## Basic Information

<b>Catalog Number:</b> 16545-1-AP	<b>GenBank Accession Number:</b> BC070123	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 220 µg/ml	<b>GeneID (NCBI):</b> 3565	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P05112	
<b>Isotype:</b> IgG	<b>Full Name:</b> interleukin 4	
<b>Immunogen Catalog Number:</b> AG9871	<b>Calculated MW:</b> 153 aa, 17 kDa	

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
ELISA, IF, WB

**Species Specificity:**  
human

**Cited Species:**  
human, rat

## Background Information

IL4 is a cytokine produced by CD4+ T cells in response to helminthes and other extracellular parasites. It promotes the proliferation and differentiation of antigen presenting cells. IL4 also plays a pivotal role in antibody isotype switching and stimulates the production of IgE. This cytokine has been applied in the treatment of autoimmune disorder like multiple myeloma, cancer, psoriasis, and arthritis. IL4 has also been extensively applied to inhibit detrimental effect of Th1. It may promote the growth of epithelial tumors by mediating increased proliferation and survival.

## Notable Publications

Author	Pubmed ID	Journal	Application
Marianne R Perera	36232315	Int J Mol Sci	IF
Luodan Yang	29056360	Exp Neurol	ELISA
Honghao Yu	35864548	J Transl Med	WB

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data