

For Research Use Only

RPL28 Polyclonal antibody, PBS Only

Catalog Number: 16649-1-PBS



Basic Information

Catalog Number:

16649-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10036

GenBank Accession Number:

BC010173

GeneID (NCBI):

6158

UNIPROT ID:

P46779

Full Name:

ribosomal protein L28

Calculated MW:

137 aa, 16 kDa

Observed MW:

16 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 16649-1-AP (RPL28 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 16649-1-PBS in a different storage buffer formulation.