

For Research Use Only

NOC4L Polyclonal antibody

Catalog Number: 17025-1-AP **1 Publications**



Basic Information

Catalog Number:

17025-1-AP

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10313

GenBank Accession Number:

BC071958

GeneID (NCBI):

79050

UNIPROT ID:

Q9BVI4

Full Name:

nucleolar complex associated 4 homolog (S. cerevisiae)

Calculated MW:

516 aa, 59 kDa

Observed MW:

59 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : HeLa cells, Jurkat cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
James C Murphy	36653366	Nat Commun	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 17025-1-AP (NOC4L antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.