

For Research Use Only

# DNAJC27 Polyclonal antibody

Catalog Number: 17043-1-AP

Featured Product

1 Publications



## Basic Information

Catalog Number:

17043-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10586

GenBank Accession Number:

BC034049

GeneID (NCBI):

51277

UNIPROT ID:

Q9NZQ0

Full Name:

DnaJ (Hsp40) homolog, subfamily C, member 27

Calculated MW:

273 aa, 31 kDa

Observed MW:

31 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human

Positive Controls:

WB : mouse brain tissue, HeLa cells, rat brain tissue

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Yujin Wang	36397632	Curr Pharm Des	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 17043-1-AP (DNAJC27 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.