

For Research Use Only

UBE2U Polyclonal antibody

Catalog Number: 17323-1-AP **1 Publications**



Basic Information

Catalog Number:

17323-1-AP

Size:

350 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11302

GenBank Accession Number:

BC029895

GeneID (NCBI):

148581

UNIPROT ID:

Q5VVX9

Full Name:

ubiquitin-conjugating enzyme E2U
(putative)

Calculated MW:

226aa, 27 kDa; 321aa, 38 kDa

Observed MW:

38 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse

Background Information

UBE2U (Ubiquitin-conjugating enzyme E2 U) is also named as Ubiquitin carrier protein U, Ubiquitin-protein ligase U and belongs to the ubiquitin-conjugating enzyme family. It catalyzes the covalent attachment of ubiquitin to other proteins. It has 2 isoforms produced by alternative splicing with the molecular weight of 38 kDa and 27 kDa.

Notable Publications

Author	Pubmed ID	Journal	Application
Yingying Guo	27903633	J Biol Chem	

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data