

For Research Use Only

TIMP2 Polyclonal antibody

Catalog Number: 17353-1-AP **39 Publications**



Basic Information

Catalog Number: 17353-1-AP	GenBank Accession Number: BC052605	Purification Method: Antigen affinity purification
Size: 700 ug/ml	GeneID (NCBI): 7077	Recommended Dilutions: IHC 1:50-1:500
Source: Rabbit	UNIPROT ID: P16035	
Isotype: IgG	Full Name: TIMP metalloproteinase inhibitor 2	
Immunogen Catalog Number: AG10051	Calculated MW: 220 aa, 24 kDa	

Applications

Tested Applications: IHC, ELISA	Positive Controls: IHC : human kidney tissue, mouse lung tissue, human prostate cancer tissue
Cited Applications: IHC, IF	
Species Specificity: human, mouse, rat	
Cited Species: human, mouse, rat, bovine	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

The TIMPs are natural highly specific endogenous inhibitors of MMPs. Human TIMP family consists of four 21-28 kDa proteins known as TIMP1, TIMP2, TIMP3, and TIMP4 encoded by four paralogous genes. TIMPs are vital to the maintenance of ECM homeostasis primarily through their MMP inhibitory functions. TIMP2, a member of the TIMP family, regulates the proteolytic activity of all MMPs and is involved in cell differentiation, growth, migration, angiogenesis, and apoptosis. Urinary TIMP2 has been recently recognized as an early biomarker to predict Acute kidney injury (AKI) in critically ill patients.

Notable Publications

Author	Pubmed ID	Journal	Application
Xiao-Xiao He	30219235	Biochem Biophys Res Commun	
Jen-Jie Lin	25222667	Mar Drugs	
Qianqian Zhang	36310322	Oncogene	

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 17353-1-AP (TIMP-2 antibody) at dilution of 1:200 (under 10x lens).




Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using 17353-1-AP (TIMP-2 antibody) at dilution of 1:200 (under 40x lens).



Immunohistochemical analysis of paraffin-embedded mouse lung tissue slide using 17353-1-AP (TIMP-2 antibody) at dilution of 1:400 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse lung tissue slide using 17353-1-AP (TIMP-2 antibody) at dilution of 1:400 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human prostate cancer tissue slide using 17353-1-AP (TIMP-2 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).