

For Research Use Only

TLX3 Polyclonal antibody

Catalog Number: 17454-1-AP



Basic Information

Catalog Number:

17454-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10612

GenBank Accession Number:

BC017291

GeneID (NCBI):

30012

UNIPROT ID:

O43711

Full Name:

T-cell leukemia homeobox 3

Calculated MW:

291 aa, 32 kDa

Observed MW:

32 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

Applications

Tested Applications:

IP, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

IP: Jurkat cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-TLX3 (IP:17454-1-AP, 4ug;
Detection:17454-1-AP 1:300) with Jurkat cells
lysate 2800ug.