

For Research Use Only

KLHL32 Polyclonal antibody

Catalog Number: 17775-1-AP



Basic Information

Catalog Number: 17775-1-AP	GenBank Accession Number: BC071627	Purification Method: Antigen affinity purification
Size: 500 µg/ml	GeneID (NCBI): 114792	Recommended Dilutions: IF/ICC 1:200-1:800
Source: Rabbit	UNIPROT ID: Q96NJ5	
Isotype: IgG	Full Name: kelch-like 32 (Drosophila)	
Immunogen Catalog Number: AG12196	Calculated MW: 470 aa, 53 kDa	

Applications

Tested Applications: IF/ICC, ELISA	Positive Controls: IF/ICC : HeLa cells,
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using KLHL32 antibody (17775-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).