

For Research Use Only

SPRY3 Polyclonal antibody

Catalog Number: 17932-1-AP



Basic Information

Catalog Number:

17932-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12364

GenBank Accession Number:

BC074973

GeneID (NCBI):

10251

UNIPROT ID:

O43610

Full Name:

sprouty homolog 3 (Drosophila)

Calculated MW:

288 aa, 31 kDa

Observed MW:

31-33 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : RAW 264.7 cells, RAW264.7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



RAW264.7 cells were subjected to SDS PAGE followed by western blot with 17932-1-AP (SPRY3 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.