

For Research Use Only

WNT3 Polyclonal antibody

Catalog Number: 17983-1-AP **4 Publications**



Basic Information

Catalog Number:

17983-1-AP

Size:

160 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12399

GenBank Accession Number:

BC114219

GeneID (NCBI):

7473

UNIPROT ID:

P56703

Full Name:

wingless-type MMTV integration site family, member 3

Calculated MW:

355 aa, 40 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human, mouse, pig

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Qing Zuo	27853262	Sci Rep	WB
Jianqing Ye	30431120	Oncol Rep	WB
Min Li	27247798	J Ophthalmol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data