

For Research Use Only

APOA5 Polyclonal antibody

Catalog Number:18019-1-AP

2 Publications



Basic Information

Catalog Number:

18019-1-AP

Size:

260 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12568

GenBank Accession Number:

BC101787

GeneID (NCBI):

116519

UNIPROT ID:

Q6Q788

Full Name:

apolipoprotein A-V

Calculated MW:

366 aa, 41 kDa

Observed MW:

42-45 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:300-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB, IHC

Species Specificity:

human, pig

Cited Species:

human

Positive Controls:

WB : HepG2 cells, pig liver tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Zhichao Lai	38454421	J Transl Med	IHC
Takashi Suzuki	38346769	Endocr J	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cell lysates were subjected to SDS PAGE followed by western blot with 18019-1-AP (APOA5 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.