

For Research Use Only

KIRREL3 Polyclonal antibody

Catalog Number: 18055-1-AP



Basic Information

Catalog Number:

18055-1-AP

Size:

133 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12573

GenBank Accession Number:

BC101801

GeneID (NCBI):

84623

UNIPROT ID:

Q8IZU9

Full Name:

kin of IRRE like 3 (Drosophila)

Calculated MW:

778 aa, 85 kDa

Observed MW:

65-70 kDa, 125 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse testis tissue, mouse brain tissue, mouse kidney tissue, rat brain tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 18055-1-AP (KIRREL3 antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.