

For Research Use Only

CENPP Polyclonal antibody

Catalog Number: 18096-1-AP



Basic Information

Catalog Number:

18096-1-AP

Size:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11865

GenBank Accession Number:

BC071726

GeneID (NCBI):

401541

UNIPROT ID:

Q6IPU0

Full Name:

centromere protein P

Calculated MW:

288 aa, 33 kDa

Observed MW:

33 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Jurkat cells, HEK-293T cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 18096-1-AP (CENPP antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.