

For Research Use Only

NFATC1 Polyclonal antibody

Catalog Number: 18337-1-AP

7 Publications



Basic Information

Catalog Number:

18337-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13014

GenBank Accession Number:

BC104753

GeneID (NCBI):

4772

UNIPROT ID:

O95644

Full Name:

nuclear factor of activated T-cells,
cytoplasmic, calcineurin-dependent 1

Calculated MW:

716aa, 78 kDa; 943aa, 101 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IF

Species Specificity:

human, mouse, rat

Cited Species:

mouse, rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Zhihui Kuang	40619563	Ann Biomed Eng	WB
Zhengchao Wang	40159624	J Cell Mol Med	IF
Kaiqiang Yang	40059399	Redox Rep	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

