

For Research Use Only

GAL3ST1 Polyclonal antibody

Catalog Number: 18503-1-AP



Basic Information

Catalog Number:

18503-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13353

GenBank Accession Number:

BC014649

GeneID (NCBI):

9514

UNIPROT ID:

Q99999

Full Name:

galactose-3-O-sulfotransferase 1

Calculated MW:

49 kDa

Observed MW:

48 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse testis tissue, HEK-293 cells, SGC-7901 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 18503-1-AP (GAL3ST1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.