

For Research Use Only

PCOTH Polyclonal antibody

Catalog Number: 18808-1-AP



Basic Information

Catalog Number:

18808-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC142999

GeneID (NCBI):

542767

UNIPROT ID:

Q58A44

Full Name:

prostate collagen triple helix

Calculated MW:

107 aa, 11 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

PCOTH may be involved in growth and survival of prostate cancer cells through the TAF-Ibeta pathway.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data