

For Research Use Only

PTRF Polyclonal antibody

Catalog Number:18892-1-AP

Featured Product

32 Publications



Basic Information

Catalog Number:

18892-1-AP

Size:

550 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13486

GenBank Accession Number:

BC008849

GeneID (NCBI):

284119

UNIPROT ID:

Q6NZI2

Full Name:

polymerase I and transcript release factor

Calculated MW:

43 kDa

Observed MW:

50-55 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:10000

IP 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

IHC 1:500-1:2000

IF/ICC 1:200-1:800

Applications

Tested Applications:

WB, IHC, IF/ICC, IP, ELISA

Cited Applications:

WB, IHC, IF, IP

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, monkey

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : HT-1080 cells, rat heart tissue, mouse heart tissue, mouse lung tissue

IP : HT-1080 cells,

IHC : human normal colon, mouse brown adipose tissue

IF/ICC : HepG2 cells,

Background Information

Polymerase I and transcript release factor(PTRF) is an essential factor in the biogenesis of caveolae, which involve in numerous process, such as signal transduction and membrane and lipid trafficking. PTRF can terminate transcription of RNA Pol I, which involves pausing of transcription by TTF1, and dissociation of the transcription complex, release the complex from the template. The calculated molecular weight of PTRF is 43 kDa, but modified PTRF is about 50-55kDa. (PMID: 11139612)

Notable Publications

Author	Pubmed ID	Journal	Application
Zhiwei Liao	34476937	ACS Nano	WB
Xiang-Li Bai	29262615	Oncotarget	WB
Xiaotong Cen	36443833	Cell Commun Signal	IP,WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 18892-1-AP (PTRF antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.




IP result of anti-PTRF (IP:18892-1-AP, 4ug; Detection:18892-1-AP 1:2000) with HT-1080 cells lysate 3600 ug.



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using PTRF antibody (18892-1-AP) at dilution of 1:400 and Coralite@488-Conjugated Goat Anti-Rabbit IgG(H+L), CL594-Phalloidin (red).



Immunohistochemical analysis of paraffin-embedded human normal colon slide using 18892-1-AP (PTRF antibody) at dilution of 1:1000 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human normal colon slide using 18892-1-AP (PTRF antibody) at dilution of 1:1000 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).