

For Research Use Only

TGF Beta 1 Polyclonal antibody

Catalog Number: 18978-1-AP **90 Publications**



Basic Information

Catalog Number: 18978-1-AP	GenBank Accession Number: NM_000660	Purification Method: Antigen affinity purification
Size: 220 µg/ml	GeneID (NCBI): 7040	Recommended Dilutions: WB 1:200-1:1000 IHC 1:50-1:500
Source: Rabbit	UNIPROT ID: P01137	
Isotype: IgG	Full Name: transforming growth factor, beta 1	
	Calculated MW: 44 kDa	
	Observed MW: 25 kDa	

Applications

Tested Applications: IHC, WB, ELISA	Positive Controls: WB : MCF-7 cells, IHC : human lung cancer tissue,
Cited Applications: WB, IHC, IF, Cell treatment	
Species Specificity: human, mouse, rat	
Cited Species: human, mouse, rat, pig, canine, monkey, sheep	
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

TGFB, also named as LAP and TGFB1, is a multifunctional peptide that controls proliferation, differentiation, and other functions in many cell types. TGFB acts synergistically with TGFA in inducing transformation. It also acts as a negative autocrine growth factor. Dysregulation of TGFB activation and signaling may result in apoptosis. Many cells synthesize TGFB and almost all of them have specific receptors for it. TGFB positively and negatively regulates many other growth factors. It plays an important role in bone remodeling as it is a potent stimulator of osteoblastic bone formation, causing chemotaxis, proliferation and differentiation in committed osteoblasts. It is highly expressed in bone. Mutation of TGFB are the cause of Camurati-Engelmann disease (CED) which known as progressive diaphyseal dysplasia 1 (DPD1). Full length, inactive 44 kD TGFB1 is cleaved into mature TGFB1 (13 kD). TGFB1 also homodimerizes and heterodimerizes with TGFB2, so there is potential for multiple different band sizes in WB (12, 25, 45 to 65 kDa).

Notable Publications

Author	Pubmed ID	Journal	Application
Fang Dou	30215298	Rejuvenation Res	WB
Xin Liang	28850927	Biomed Pharmacother	WB
Haoyu Ruan	27581744	Sci Rep	IHC

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



MCF-7 cells were subjected to SDS PAGE followed by western blot with 18978-1-AP (TGF-beta 1 antibody at dilution of 1:300 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using 18978-1-AP (TGF-beta 1 antibody at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using 18978-1-AP (TGF-beta 1 antibody at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).