

For Research Use Only

TLR1-Specific Polyclonal antibody

Catalog Number: 19500-1-AP **1 Publications**



Basic Information

Catalog Number: 19500-1-AP	GenBank Accession Number: NM_003263	Purification Method: Antigen affinity purification
Size: 130 µg/ml	GeneID (NCBI): 7096	
Source: Rabbit	UNIPROT ID: Q15399	
Isotype: IgG	Full Name: toll-like receptor 1	
	Calculated MW: 90 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
WB

Species Specificity:
human

Cited Species:
rat

Background Information

TLR1, also named as TIL, CD281 and KIAA0012, belongs to the Toll-like receptor family. TLR1 participates in the innate immune response to microbial agents. TLR1 cooperates with TLR2 mediates the innate immune response to bacterial lipoproteins or lipopeptides. Acts via MYD88 and TRAF6, leading to NF-kappa-B activation, cytokine secretion and the inflammatory response. This antibody is specific to TLR1.

Notable Publications

Author	Pubmed ID	Journal	Application
Dongmin Li	25534427	Mol Cell Endocrinol	WB

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data