

For Research Use Only

# SPACA5 Polyclonal antibody

Catalog Number: 19525-1-AP



## Basic Information

Catalog Number:

19525-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_205856

GeneID (NCBI):

389852

UNIPROT ID:

Q96QH8

Full Name:

sperm acrosome associated 5

Calculated MW:

18 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

SPACA5, also named as LYZL5, SPACA5A and SPACA5B, belongs to the glycosyl hydrolase 22 family. SPACA5 hydrolysis (1→4)-beta-linkages between N-acetylmuramic acid and N-acetyl-D-glucosamine residues in a peptidoglycan and between N-acetyl-D-glucosamine residues in chitodextrins.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

