

For Research Use Only

UVRAG Polyclonal antibody

Catalog Number: 19571-1-AP **4 Publications**



Basic Information

Catalog Number:

19571-1-AP

Size:

260 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13765

GenBank Accession Number:

BC064837

GeneID (NCBI):

7405

UNIPROT ID:

Q9P2Y5

Full Name:

UV radiation resistance associated

gene

Calculated MW:

699 aa, 78 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB, IP

Species Specificity:

human, mouse, rat

Cited Species:

human, rat

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Cefan Zhou	31696776	Autophagy	WB
Vinod Kumar	31116585	Am J Physiol Cell Physiol	WB,IP
Hisao Naito	31467308	Sci Rep	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data