

For Research Use Only

PLP2 Polyclonal antibody

Catalog Number: 19710-1-AP



Basic Information

Catalog Number: 19710-1-AP	GenBank Accession Number: NM_002668	Purification Method: Antigen affinity purification
Size: N/A	GeneID (NCBI): 5355	
Source: Rabbit	UNIPROT ID: Q04941	
Isotype: IgG	Full Name: proteolipid protein 2 (colonic epithelium-enriched)	
	Calculated MW: 17 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

PLP2, also named as A4, may play a role in cell differentiation in the intestinal epithelium.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data