

For Research Use Only

ZNF26 Polyclonal antibody

Catalog Number: 19839-1-AP



Basic Information

Catalog Number:

19839-1-AP

Size:

550 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13837

GenBank Accession Number:

BC008211

GeneID (NCBI):

7574

UNIPROT ID:

P17031

Full Name:

zinc finger protein 26

Calculated MW:

61 kDa

Observed MW:

70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, HepG2 cells, MCF-7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 19839-1-AP (ZNF26 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.