

For Research Use Only

REV3L Polyclonal antibody

Catalog Number: 19985-1-AP



Basic Information

Catalog Number:

19985-1-AP

Size:

N/A

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002912

GeneID (NCBI):

5980

UNIPROT ID:

O60673

Full Name:

REV3-like, catalytic subunit of DNA polymerase zeta (yeast)

Calculated MW:

353 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human, mouse, rat

Background Information

REV3L, also named as POLZ and REV3, belongs to the DNA polymerase type-B family. REV3L catalyses the reaction: Deoxynucleoside triphosphate + DNA(n) = diphosphate + DNA(n+1). The antibody recognizes the C-term of REV3L.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

